




The Open

The Orchid Spec Open-sided range for hysteroscopic, sonographic, colposcopic and ablation procedures allows for optimal lateral access. Increase overall access and optionally leave instruments or tubing in place while removing the speculum for patient comfort.

Open-sided speculum

Allows for easy speculum retraction while leaving instruments or catheters engaged with the cervix.
Provides easy vaginal wall access.



Ideally suited for:

- IUD insertion
- Colposcopy
- Cervical biopsy
- Sonohysterography (SIS / GIS)
- HyFoSy / HyCoSy
- Office hysteroscopy
- HSG / THL
- Endometrial ablation
- LEEP / LLETZ (Open SPX model)
- Embryo transfer
- Sidewall access

Robust design that does not cave


The robust hinge and bill design allows for a spacious and stable side opening even under significant load.



reddot design award
winner 2013 & 2016

A better speculum, a better practice, a better world

The Open



Slim organic design

Its friendly, 'Clean White' appearance puts the patient more at ease, reducing pelvic tension. The slim, organic, ergonomic design makes insertion effortless, resulting in less anxiety and faster and easier procedures.

A better speculum

1. Anti-pinching gap

Due to a large gap, the bills do not touch when closing.

2. Less dripping of discharge

Due to the inward folded edges that retain fluids or discharge. Less cleaning, more hygienic, quicker procedure times.

3. Backward tilted handle

The 105° backward angled handle allows for more freedom of movement, greater range of access, deeper insertion, less changing to other speculum sizes and reduces unintended patient contact.

4. Soft rounded edges

There is no longer a risk of tissue trauma or bleeding due to sharp parting or product lines. Our patented design enables uniquely rounded, smooth edges both on the in and outside edge, improving patient comfort significantly.

5. Reflective surface

The white surface reflects and directs external light to the tip much better than a metal or transparent speculum. An inner light source is not needed.

6. 30% wider opening

compared to most used specula. Provides a better view and allows for easier manipulation of instruments.

7. Click-free operation

Easy single-handed (un)locking without any fidgeting, rattling or other noises that cause unwanted tensing of the patient.

8. Easy single-handed use & no risk of jamming

Designed for single-handed operation, including easy locking and unlocking. The lock has a safety mechanism that will close the bills if too much force is applied.

“The Orchid Spec specula are in a class of their own; my patients say their examination is more comfortable. The Open SPX is my treatment speculum of choice with improved access and visual field.”

Prof. Maggie Cruickshank,
President of the European Federation for Colposcopy
Consultant Gynaecologist
Emeritus Professor University of Aberdeen



The Open SPX

Smoke plume extraction (optional)




The tapered nozzle allows for easy connection of all standard smoke suction tubes and provides a reliable connection.



Enhanced cervix support

The patented inward folded edges combined with the flat support surfaces at the tip of the bills allow for improved lifting and positioning of the cervix.

We offer the Open-sided in two sizes: Medium and Wide. Both sizes come in our popular 'Clean White' and 'Crystal Clear'. The smoke plume extraction version is available in 'Clean White' Medium and -Wide.



Use a white or clear speculum?

Our White speculum

Reduced misdiagnosis

Diagnosing through a 'clear' speculum will provide a distorted view due to local tissue depression and resulting tissue discoloration.

Improved examination

The white surface permits quicker and better examination of discharge color.

Patient friendlier


Patients perceive white as higher quality, hygienic and a friendlier appeal. Less anxious patients improve procedure times and efficiency.

Preferred by 95% of our clients

even though we also sell a clear speculum.

Our Clear speculum

The Orchid Spec's clear speculum provides a visible advantage in accessing larger conditions, bleeding, or inspecting sutures, for instance, after a hysterectomy.



"The best speculum available today. Both patient- and user-friendly, with optimal performance during all procedures indicated."

Prof. Dr. M. H. Emanuel,
Gynecologist at UMC Utrecht, and visiting professor at the university hospital of Ghent.


Today, the Bio is
100%
sugarcane

100%
carbon neutral
by 2025

net
negative
by 2030

At Orchid Spec, we aspire to make vaginal procedures as effortless as possible for women and gynecologists worldwide - while helping mother nature. The Orchid Spec is the most ecological choice today. We aim to be CO₂ neutral by 2025 and Net Negative by 2030.


All our products are Dutch Design and manufactured in The Netherlands in a state-of-the-art facility to guarantee optimal management of natural resources and waste reduction.



The Standard
Our standard
ecological speculum



The Bio
Made from
100% sugarcane



The Open
With side-opening



The SPX
With surgical
plume extraction



The Open SPX
With side-opening and
surgical plume extraction

A better speculum, a better practice, a better world



For more information or enquiries, please contact us at:

www.orchid-spec.com

enquiry@brideamedical.com

The Orchid Spec products are patented and design registered. Orchid Spec and Bridea Medical are registered trademarks of Bridea Medical from The Netherlands. © 2025 Bridea Medical.